

Molecular Devices SpectraMax M5 Multi-Mode Microplate Reader
Format Documentation for SoftMax Pro Export Files.

File Format:

- Export starts with ##BLOCKS= followed by number of data blocks.
- Each block starts with Plate: header containing tab-separated metadata.
- Metadata fields: plate name, version, PlateFormat, read mode (Endpoint/Kinetic), detection method (Fluorescence/Absorbance/Luminescence), data type (Raw), autoscale flag, excitation wavelength (nm), emission wavelength (nm), first column, last column, number of wells, first row, last row.
- Data grid: Temperature row followed by rows A through H.
- Each data row: row letter, temperature in Celsius, then values for columns 1-12.
- Block ends with ~End marker.

Detection Methods:

- Fluorescence Intensity: requires excitation and emission wavelengths.
- Absorbance: single measurement wavelength, no emission.
- Luminescence: no excitation or emission required.

Plate Formats: 96-well, 384-well.

Well naming: row letter (A-H for 96-well) + column number (1-12 for 96-well).

Values are raw fluorescence units (RFU) for fluorescence mode.

Temperature in degrees Celsius measured at read time.